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

DIGITAL pH METER

PORTABLE, TABLE STAND, WALL MOUNT

Before putting the pH meter to use, please read the following instructions carefully.

- (2) Turn on the power switch of instrument.
- (3) Waiting few seconds for stable reading, adjust the "CAL.(pH7) adj. knob" until the display reading same as 7.00 exactly.
- (4) Rinse the electrode in distilled water and drying it.
- (5) Place electrode into buffer pH4.00 solution. Waiting few seconds for stable reading, adjust "SLOPE (pH4) adj. knob" until the display reading same as 4.00 exactly.
- (6) Repeat above procedures two times at least.

LCD BACK LIGHT

It will automatically turn on the LCD back light when plug in the AC adaptor power supply. (Attention the DC power polarity).

MEASURING PROCEDURE


After the instrument and pH electrode are calibrated, then the unit is now ready for measuring.

- (1) Connect the combination pH electrode to the BNC socket.
- (2) Turn on the power switch of instrument.
- (3) Place the pH electrode into the measured solution, then the instrument will display the pH value.
- (4) After make the measurement, please rinse the electrode in distilled water. Then put the electrode into bottle with 3M KCL buffer solution for storage.

REPLACEMENT OF BATTERY

- (1) When the left corner of LCD display battery empty, it indicate a normal battery output of less than 6.5 - 7.5V. It is necessary to replace the battery.
- (2) Replace the battery, remove the battery cover (back).
- (3) Take out the battery, install a new one (DC9V 006P, 6F22K) and reinstate the battery cover again.

FRONT PANEL DESCRIPTION

- | | |
|----------------------------|----------------------------------|
| 1) LCD Display | 6) pH electrode BNC input socket |
| 2) Low battery indicator | 7) Stand |
| 3) Power ON / OFF switch | 8) Battery cover |
| 4) CAL. (pH7) adjust knob | 9) AC regulator adapter socket |
| 5) SLOPE (pH4) adjust knob | (2.5 mm shaft) |

PH CALIBRATING PROCEDURE

Calibrating Consideration

These pH meter already calibrated by mV signal that simulated from the ideal pH Electrode mV output (base on 25 °C environment). However due to (1) An ideal electrode will produce 0 mV at pH7.00, but most electrodes are slightly off. (2) The measuring environment temperature may not near 25 °C. It is necessary to make the following calibration procedure that to keep instrument combine electrode within high accuracy.

Requiring Equipment for Calibration

- (1) Combination pH Electrode (optional).
- (2) Two buffer solutions (optional) : pH7.00 & pH4.00.

Two Points Calibration

- (1) Connect the pH Electrode to the BNC socket, rinse it in distilled water and drying it, place electrode into buffer pH7.00 solution. (Avoid pollutes the calibrating buffer solution.)

INSTRUCTION FOR CLEANING AND MAINTENANCE

- 1) It must to turn off the power when cleaning the instrument.
- 2) Please wipe the housing with drying when cleaning, do not wipe with wetting.
- 3) It had better calibrate the pH electrode during 2-4 weeks to make sure the reading within accuracy when measuring in a long time.

GENERAL SPECIFICATION

Display : 15mm (0.59") LCD 3 1/2 digits.
 Measuring Range : 0 to 14 pH.
 Resolution : 0.01 pH.
 Accuracy : $\pm(0.1\%+2 \text{ digits})$ after calibration.
 Input Impedance : 10^{12} ohms.
 Sampling Time : 0.4 seconds.
 Operating Temp. : 0 to 50 °C. (32 to 122 °F).
 Operating Humidity : Less than 90% RH.
 Calibration Knob : External, pH4 (Slope), pH7 (Cal.)
 Power Supply : DC9V 006P, 6F22K battery,
 AC power adapter (optional)
 Power Consumption : Approx. DC 2 mA. (LCD back light off).
 Dimension : 110 x 60 x 25 mm. (4.3 x 2.4 x 1 inch).
 Weight : Approx. 98 g. (0.2 lb.) without battery.
 Standard Accessories : Instruction Manual, Adjustment Sticker,
 Electrode Holder, pH Meter Mounting Holder.
 pH Electrode (optional) : Any pH Electrode with BNC connector.
 AC regulator adapter (optional) : output DC9V 125 mA.
 (2.5mm shaft)